




Zirc Dental Products

Product Catalog

Helping dental teams be safe, happy and efficient.

2 - 12 COLOR METHOD

- 2 Build Your Own Tub
- 3 All-In-One Kits
- 4 Tub Accessories
- 5 Trays & Covers
- 6 Multi-Mod Racks
- 7 Cassettes & Instrument Holders
- 8 Performance Tool Organizers & Steri-Containers
- 9 Bur Blocks & Covers
- 10 Drawer Organizers, Steri-Bur Guards & Adapters
- 11 Steri-Endo Guards & Holders
- 12 Rings & Tape

13 - 14 ISOLATION

- 13 Airway Armor®, Mr. Thirsty®, and SE Cushion
- 14 Insti-Dam®, Pink Petal™, and Lingua-Fix®

15 - 17 MOUTH MIRRORS

- 16 Crystal HD® Mirrors
- 17 Mirror Handles & Mirror Magic® Anti-Fog Solution

18 - 23 OPERATORY/MISC

- 18 Denture & Retainer Boxes
- 19 Steri-Soaker, CAD/CAM Block Organizer, and Instrument Mats
- 20 E-Z Storage & Handi-Hopper
- 21 E-Z Shelf & Accessories
- 22 Water Wise & Dispensers
- 23 Mighty Mixer Bowls, Spatula, and Organizers





colormethod TRANSFORM EFFICIENT INTO PROFIT

Maximize your team's potential and watch profits soar with instrument and material management.

Benefits include:

- ✓ **Efficient Teams**
- ✓ **Higher Profits**
- ✓ **Lower Costs**
- ✓ **Optimized Workspaces**
- ✓ **Smoother Onboarding**



Experts that can help

Learn more at ColorMethod.com or call (763) 682-6636

Scan to see what people say about us.



Build Your Own Tub



- 1** 20Z463
Procedure Tub
\$23 Available in all colors
Allows you to "build" your own setup.
Dimensions: 12-1/4" x 10-7/8" x 2-3/4"
S

- 2** 20Z462
Tub Liner
\$6.50
Foam liner that sits in bottom of Tub.
Tub not included. Replace as needed.
Dimensions: 9-5/8" x 8-7/8" x 1/8"

- 3** 20Z461
Tub Dividers
\$15.50 Gray only
Creates up to 24 compartments. May be cut to create larger compartments.
Dimensions: 9-5/8" x 8-5/8" x 1-5/16"
S

- 4** 20Z459
Safe-Lok® Tub Cover
\$23
Locking cover for securing materials.
Keeps materials protected during storage and transportation.
Dimensions: 12-5-8" x 11-1/4" x 1-7/16"
S

- 5** 20Z465
Tub Insert
\$18 Gray only
10 compartments for easy organization. Sits on rim of Tub, allowing for storage underneath.
Dimensions: 10-3/4" x 9-11/16" x 1"
S

- 6** 20Z474
Tub Box
\$16.50 Gray only
Designed to hold/organize supplies.
Non-slip bottom for secure placement.
Flip-top cover.
Dimensions: 3-3/4" x 3-5/8" x 1-7/8"
S

All-In-One Kits

Zirc's procedure tubs are designed to store and organize materials to keep your inventory under control.



20Z455

Complete Tub

\$87 Available in all colors

Includes:

- Procedure Tub
- Tub Liner
- Tub Dividers
- Single Tub Cup/Cover
- Double Tub Cup/Cover
- Divided Slide Tray
- Safe-Lok® Cover for Tub

Dimensions: 12-5/8" x 11-1/4" x 4"

S

Composite Kits S

Store and organize composite materials with our convenient Composite Kits.



20Z489

Syringe Composite Kit

\$76

Syringe holder allows you to store up to 40 syringe tubes.

Includes:

- Syringe 20-Unit Stand
- 1 Long Tub Cup/Cover
- 1 Divided Slide Tray
- 1 Syringe Clip-On Holder
- Shade Label Sheet 2pk



20Z479


Capsule Composite Kit

\$64

Long Tub Cups provide a total of 20 compartments for capsule storage.

Includes:

- Tub Liner
- 4 Long Tub Cup/Cover
- 1 Divided Slide Tray
- 1 Capsule Clip-On Holder
- Shade Label Sheet 2pk



Tub Accessories

Tub Cups with Covers S

Fill Tub Cups with brushes, clamps, cotton pellets, prophylaxis angle tips, prophylaxis paste, pumice, syringe tips, implant tools, etc. Designed to be used inside the Procedure Tub.



20Z473
Long Tub Cup w/ Cover
\$16 Gray only

Internal dividers are removable.

Dimensions: 8-3/8" x 2-5/16" x 1-1/8"



20Z472
Double Tub Cup w/ Cover
\$10 Gray only

Dimensions: 3-3/8" x 2-1/4" x 1-1/8"



20Z471
Single Tub Cup w/ Cover
\$10 Gray only

Dimensions: 2-1/4" x 1-5/8" x 1-1/8"




20Z477
Shade Label Sheets 2pk
\$6

Preprinted labels for marking composite materials within the Tub organizers.

Dimensions: 7" x 4-3/8" x 1/64"


Syringe Holders S W



20Z487
Syringe 20-Unit Stand
\$26

Stands up for easy access to syringe-able materials. Lays flat to fit inside Tub or drawer. Holds up to 20 syringes. Has space for 2 Syringe 10-Unit Add-Ons (Up to 40 syringes).

Dimensions: 8-3/4" x 4-3/4" x 1-9/16"



20Z486
Syringe Clip-On Holder
\$11

Allows quick access to composite during procedure. Conveniently clips on any side of the Tub or B-Lok Tray. Holds 15 syringes and 3 capsules.

Dimensions: 8" x 2-7/16" x 3/4"



20Z488
Syringe 10-Unit Add-On
\$11

An add-on for the Syringe 20-Unit Stand. Add for storage of up to 40 syringe tubes. Gives additional space for syringe storage. Holds 10 syringes.

Dimensions: 8-9/16" x 2-5/16" x 3/4"



20Z476
Capsule Clip-On Holder
\$11

Allows quick access to composite during procedure. Conveniently clips on any side of the Tub or B-Lok Tray. Holds 8 capsules, and 4 syringes.

Dimensions: 8" x 2" x 9/16"

Slide Trays S W

Sits on rim of Tub and allows for storage underneath.



20Z466
Divided Slide Tray
\$15 Clear

Dimensions: 9-3/4" x 5-3/8" x 1-1/16"



20Z467
Flat Slide Tray
\$15 Clear

Dimensions: 9-3/4" x 5-3/8" x 1-1/16"

Trays & Covers

Our exclusive B-Lok trays and covers provide an excellent method of infection control and safety that allows for two benefits other trays do not: preventing any contaminated items from falling off the used tray during transportation from operatory to central sterilization, and keeping contaminated items enclosed until you are ready for the cleaning and sterilization process.



20Z401
B-Lok Flat Tray
\$15.50 Available in all colors

Wide edges for secure grip. Non-slip feet for stability. High-gloss surface helps prevent staining.


Dimensions: 13-3/8" x 9-5/8" x 7/8"



20Z451
B-Lok Divided Tray
\$17.50 Available in all colors

Holds 12 hand instruments, including 9 FG, RA, or HP bur holders. Features Insti-Grip, non-slip feet for stability, and self-lubricating, stain-resistant surface. Also includes holder for Syringe Tip or Saliva Ejector.

Dimensions: 13-3/8" x 9-5/8" x 7/8"



20Z454I
Insti-Grip
\$4

To be used with the B-Lok divided Tray only. Holds instruments quietly and safely in place.

Dimensions: 6-5/8" x 1/4" x 3/8"



20Z101
Mini Tray
\$14 Available in all colors

Wide edges for secure grip. Non-slip feet for stability. High-gloss surface helps prevent staining. Mini Trays are great for hygiene/exam setup.

Dimensions: 9-3/8" x 6-3/8" x 7/8"



20Z445
Safe-Lok® Tray Cover
B-Size
\$20

Fits B-Lok Flat Tray and Divided Trays only. Slide locking mechanism securely locks cover to tray.

Dimensions: 13-7/8" x 9-7/8" x 3/4"



Tall
B-Size

Standard
B-Size

20Z446
Safe-Lok® Tray Cover
Tall B-Size
\$21

The Tall size accommodates our Performance Tool Organizer and almost any stainless steel cassette. Fits B-Lok Flat and Divided Trays only. Side locking mechanism securely locks cover to tray.

Dimensions: 13-7/8" x 9-7/8" x 1-3/8"



20Z441
Tray Cover B-Size
(Non-Locking)
\$19.50

Dimensions: 14" x 10-1/8" x 5/16"



20Z103
Safe-Lok® Mini Tray Cover
F-Size
\$18

Side locking mechanism securely locks cover to tray.

Dimensions: 9-7/8" x 6-5/8" x 1"



20Z102
Mini Tray Cover F-Size
(Non-Locking)
\$17.50

Dimensions: 9-5/8" x 6-5/8" x 7/8"



Multi-Mod Racks

Tub and tray racks are effective additions to a dental practice. These racks have non-slip feet so they are stable on any counter or cabinet surface.



21Z410
Multi-Mod Rack 8-Place
\$158

Holds 8 Trays or 4 Tubs with covers. Sturdy metal construction; chrome-plated. Non-slip feet. Adjustable to hold trays and tubs.

S



21Z105
Multi-Mod Rack 6-Place
\$151

Holds 6 trays or 3 tubs with covers. Adjustable to hold trays and tubs. Sturdy metal construction; chrome-plated. Non-slip feet.

S



21Z101
Mini Tray Rack
\$120

Holds 6 mini trays with covers. Chrome-plated, sturdy metal construction with non-slip feet. Fits size F trays.

S



Are you in compliance?

The Safe-Lok® Tray Cover helps in compliance with OSHA and CDC.

OSHA

Immediately or as soon as possible after use, contaminated reusable sharp shall be placed in a puncture-resistant, leak-proof and labeled or color-coded containers until properly reprocessed. To avoid spillage of contents, it is suggested that they be covered and secured prior to moving.

CDC


Contaminated instruments must be "contained" prior to transporting them. Containment can be defined as instrument cassettes or cages in which dental care items are cleaned, sterilized, and stored until point of use. Sharp instruments should not be carried openly to the instrument-processing area. Percutaneous injuries may occur to other dental healthcare personnel (DHCP) or patients in hallways leading to the instrument-processing area.

ZIRC TIP

Instrument Holders




Individual instrument handling takes time away from treating patients and increases overall reprocessing time. Zirc's cassettes improve instrument efficiency, safety, and compliance without breaking the bank.



50Z935
E-Z Jett® Pro
\$90 Available in all colors

The E-Z Jett® Pro Cassette features a spacious, multi-layout design that can hold up to 20 hand instruments including larger specialty tools. The cassette's removable components offer flexibility, allowing for easy customization based on the specific needs of your practice. The cover can be attached and removed using the three latches on each side.


Dimensions: 9.7882" x 8.2108" x 1.5770"



50Z930
E-Z Jett® Cassette
10-Place
\$51 Available in all colors

Holds 10 instruments. Accessory compartment in the middle for extra storage. Covers stay attached to the cassette and fold back easily.

Dimensions: 7-15/16" x 7" x 1-1/8"
Sterilization Pouch Size: 10" x 13" or 14"



50Z911
E-Z Jett® Cassette
Slim 8-Place
\$41 Available in all colors

Holds 8 instruments. Covers stay attached to the cassette and fold back easily. Compact size. Fits the Statim. A/W Syringe Tip Holder included.

Dimensions: 6-7/8" x 3-9/16" x 1-3/16"
Sterilization Pouch Size: 6" x 12"


* Size comparison with the 8-Place E-Z Jett® Cassette



50Z927
E-Z Jett® Cassette
8-Place
\$46 Available in all colors

Holds 8 instruments. Covers stay attached to the cassette and fold back easily. A/W Syringe Tip Holder included.

Dimensions: 7-15/16" x 4-7/16" x 1-1/8"
Sterilization Pouch Size: 7 1/2" x 13"




50Z925
E-Z Jett® Cassette
5-Place
\$36 Available in all colors

Holds 5 instruments. Great for basic exam setup! A/W Syringe Tip Holder included. Covers stay attached to the cassette and fold back easily. Fits Statim Sterilizer.

Dimensions: 7-5/16" x 3-1/8" x 1-3/16"
Sterilization Pouch Size: 6" x 12"

50Z926
● **Syringe Tip Holder 5pk**
\$5.25

Dimensions: 1" x 3/8" x 3/8"



Performance Tool Organizer



A deeper, roomier cassette made for large and specialty instruments. Designed for superior water flow for superb cleaning. Compact styling to fit in most sterilizers. Covers stay attached to the cassette and fold back easily. Comes with two strips of 10 customizable, removable feet to design your own layout.

50Z940

Performance Tool Organizer

\$65 Available in all colors

Dimensions: 8-3/16" x 4-1/4" x 1-7/32"
Sterilization Pouch Size: 7 1/2" x 13"

50Z944

Performance Tool Organizer Replacement Feet

\$8.50 Gray only | Two strips of 10.



50Z900

Steri-Container

\$14 Available in all colors

Keeps instrument setups together through sterilization processes. Holds 10-14 instruments.

Dimensions: 8-1/8" x 1-7/8" x 1-7/8"



50Z905

Compact Steri-Container


\$14 Available in all colors

Easy open, push down latch. Keeps instrument setups together through sterilization processes. Statim 2000 holds 7 Compact Steri-Containers. Holds 8-12 instruments.

Dimensions: 7-1/8" x 1-1/2" x 1-1/2"



ZIRC TIP



Closed

Hold your large & specialty tools!

Performance Tool Organizer fits on Zirc Trays with the Safe-Lok® Tall Cover!

- Ideal for:
- Surgical elevators and forceps
 - Orthodontic pliers
 - Hemostats
 - Scissors
 - Any instrument that does not fit in your regular cassette



Bur Blocks

Our Bur Blocks will keep your burs safe and organize them during the sterilization process.

Magnetic Bur Blocks

- Fits FG, RA, and HP Burs
- Magnetic base holds burs securely
- Non-slip feet for stability
- Drain holes help prevent rusting of burs



50Z400
**Magnetic Bur Block
8-Hole**
\$16.50 Available in all colors
Covers sold separately.
Dimensions: 1-1/2" x 1-1/2" x 1/2"



50Z403
**Magnetic Bur Block
14-Hole**
\$17 Available in all colors
Covers sold separately.
Dimensions: 2-15/16" x 1-1/2" x 1/2"



Recessed area to hold used burs.



Bur Block Covers

- Covers available to maintain cleanliness of burs during storage
- 1" covers for FG and RA burs; 2" cover (14-hole only) for lab and surgical burs
- Covers are NOT autoclavable
- Surface disinfect only




50Z401
**Magnetic Bur Block
8-Hole Cover
1" High Standard**
\$8.50
Dimensions: 1-5/8" x 1-5/8" x 1-1/16"



50Z405
**Magnetic Bur Block
14-Hole Cover
2" High HP Lab Burs**
\$8.50
Dimensions: 3-1/8" x 1-5/8" x 2-1/4"



50Z404
**Magnetic Bur Block
14-Hole Cover
1" High Standard**
\$8.50
Dimensions: 3-1/8" x 1-5/8" x 1-1/16"



Drawer Organizers S

20Z210

Drawer Organizer

\$38 Color Code: A, G, I, N, P, R, S

Each Drawer Organizer will come with 1 large divider, 2 medium dividers and 4 small dividers. Additional dividers can be purchased.

Dimensions: 11.5" x 6.6" x 2"

Large Divider (Included)

Dimensions: 0.95" x 1.97" x 11.35"

Drawer Organizer dividers

\$10 Color Code: A, G, I, N, P, R, S


Dividers may be ordered separate. Size (S) comes in 2 pieces and size (M) in 1 piece.

20Z207 | Small

Dimensions: 0.63" x 1.96" x 3.19"

20Z208 | Medium

Dimensions: 0.63" x 1.97" x 6.46"



Steri-Bur Guards & Adapters

This gadget has it all! Not only does it protect your burs during cleaning and sterilization, but it can also be positioned at three different angles to give you easy access at chairside.

Steri-Bur Guard 136 Washer Disinfect

- Fits FG, RA, and HP Burs
- Adjustable height for burs of different lengths (12-Hole Tall Bur Guard not adjustable)
- Holds burs in position for maximum cleaning effectiveness
- Drain holes help prevent rusting of burs



50Z410 Steri-Bur Guard 22-Hole

\$18 Available in all colors

Accommodates burs from 12mm-30mm.

Dimensions: 4-7/8" x 7/16" x 1-1/2"



50Z408 Steri-Bur Guard Tall 12-Hole

\$17.50 Gray only

Accommodates burs from 30.9mm-37.1mm.

Dimensions: 2-7/8" x 3/8" x 1-5/8"



50Z412

Universal Bur Adapter 5pk

\$5.25

Can be used with 12 or 22-Hole Bur Guard. Enables long and short burs to be used in the same bur guard without short burs falling out.

A short bur adapter is included to stop any shorter burs from falling out.

50Z406 Steri-Bur Guard 12-Hole

\$17 Available in all colors

Accommodates burs from 12mm-30mm.

Dimensions: 2-7/8" x 3/8" x 1-3/8"



Steri-Endo Guards & Holders

A unique delivery system to organize, store, and sterilize endodontic files. The Endo Guard swivels 360° with stops at every 45° for operator's convenience. Large slotting provides complete visibility and draining during cleaning.

Steri-Endo Guard

- Holds most lengths for staging, cleaning, and sterilizing files and reamers
- Millimeter gauge measures and positions rubber stopper
- Counter on locking bar tracks the number of sterilizations
- Right or left-hand use; optional angling
- Available in hand or engine styles
- Single-use foam insert cleans/holds files (packaged with 12)



50Z450 Steri-Endo Guard Hand

\$35.50 Available in all colors

Organizes up to 16 hand files/reamers.

Dimensions: 5-3/8" x 3/8" x 2-3/16"



50Z451

- **Foam Inserts 48pk**
\$35.50 Single-Use Only

50Z451V

- **Foam Inserts 144pk**
\$70 Single-Use Only

50Z452

Steri-Endo Guard Engine

\$35.50 Available in all colors

Organizes up to 11 engine and 5 hand files/reamers.

Dimensions: 5-3/8" x 3/8" x 2-3/16"



50Z460 Endo Assist

\$31.50 Available in all colors

The Endo Assist features a built-in millimeter gauge for measuring and positioning of rubber stoppers. Designed for both right and left-handed use, it includes small and large rings to fit most finger sizes. Its flat bottom allows it to stand independently, and it comes with a single-use foam insert for file cleaning and storage (packaged w/ 12).

Dimensions: 1-1/2" x 1-3/4" x 1-5/8"



50Z462 Endo Assist Stand

\$18 Black only

Holds endo files and reamers in upright position. Foam provides fast and easy cleaning of the files during use.

Dimensions: 1-7/8" x 1-7/8" x 1-1/4"



50Z466 Endo Assist Stand Foam 48pk

\$34.50 Single-Use Only

Dimensions: 10-1/2" x 7-3/4" x 7/8"

50Z466V Endo Assist Stand Foam 144pk

\$70 Single-Use Only

Dimensions: 10-1/2" x 7-3/4" x 7/8"

Rings & Tape

How to effectively identify your instrument setup

1 Color for Procedure

Vibrant Orange is for identifying the procedure



Did you know?

The diagonal setup of the orange rings shows the order of how the instruments should be used. An optional second color may be used to identify the practitioner.



ZIRC TIP

Place your instrument to determine which ring size you need.

SM

3.17mm-6.35mm

LG


6.35mm-9.52mm

XL

9.52mm-12.7mm

E-Z ID Rings

- Prepackaged in an easy-to-use dispenser
- Available in Small, Large, and Extra Large sizes to fit all instrument handles



Single Color 25pk

Available in all colors

SM \$8 | 70Z100

LG \$8.25 | 70Z200

XL \$8.50 | 70Z000

Vibrant System 200pk

Assorted colors
[N, O, P, R, Q, S]

SM \$53 | 70Z110

LG \$54 | 70Z210

XL \$55 | 70Z010

Classic System 200pk

Assorted colors
[A, B, D, G, I]

SM \$53 | 70Z105

LG \$54 | 70Z205

XL \$55 | 70Z005

Jewel System 200pk

Assorted colors
[E, J, T, I, M, U]

SM \$53 | 70Z115

LG \$54 | 70Z215

XL \$55 | 70Z015

E-Z ID Tape

- Tape dispenser has built-in tape cutter
- Dispenser roller prevents tape from twisting



E-Z ID Tape System

\$51 Includes 8 - 3' rolls

Classic | 70Z305

Assorted colors
[A, B, D, G, I]

Jewel | 70Z315

Assorted colors
[E, J, T, I, M, U]

Vibrant | 70Z310

Assorted colors
[N, O, P, Q, R, S]

70Z300

E-Z ID Tape Roll

\$12 Available in all colors

Includes 1 - 10' roll.





- Aspiration and ingestion prevention
- Patient protection and comfort
- Easy placement and removal

Airway Armor® 25pk \$25

50Z960SM | Small

50Z960MD | Medium

50Z960LG | Large

- Comes in 3 convenient sizes
- Single-Use Only
- Latex-free



50Z960TRI 15pk Trial Kit

\$20 Includes 5 of each size

(5) Small

(5) Medium

(5) Large



Protection & Easy Placement

mr. thirsty®

- Keeps work environments dry
- Fully assembled
- Latex-free

50pk \$128

● 50Z989SM | Small

● 50Z989MD | Medium

● 50Z989LG | Large

100pk \$230

● 50Z989SM-B | Small

● 50Z989MD-B | Medium

● 50Z989LG-B | Large

- Trimmable & disposable
- Affordable & Single-Use Only
- Adaptable & comfortable

50Z988 Variety Pack

\$86 Includes 10 of each size

● (10) Small

● (10) Medium

● (10) Large



Smaller size is perfect for smaller mouths!



50Z992 Mr. Thirsty® Comfort Kit \$70

The Comfort Kit's lightweight tubing allows your heavy HVE Valve/Tubing to stay mounted to the equipment.

Includes:

Lightweight lever HVE Valve, 6' lightweight tubing.


Mr. Thirsty® NOT included.

SE Cushion

Prevents tissue grabbing and irritation and is designed to quiet the noise of suction to increase patient comfort.

- Used with a Saliva Ejector
- Single-Use Only

50Z970 SE Cushion 100pk \$25



INSTI-DAM®

Insti-dam® features a built-in flexible medium gauge material that comes pre-assembled with a frame and pre-punched holes, ensuring easy positioning while keeping the eyes and nose open and comfortable.

- Additional holes can be punched
- 50% faster placement than traditional dam
- Single-Use Only
- Latex-free

50Z457

**Latex-Free Insti-Dam®
Relaxed Fit 20pk**
\$77

Ideal for posterior isolation. The extra "relaxed" material makes for quick posterior placement – no stretching, tearing, or pulling.

Dimensions: 4-7/8" x 4-7/8" x 6-1/8" diameter

Patent #7044738

50Z471

Insti-Dam® Dispenser
\$22 Only available in white

Can be mounted to wall or cabinet or can sit on a counter. Non-slip rubber bottom for stability. Holds 35 Insti-Dams®. Made for Latex-Free Insti-Dam® ONLY.

Dimensions: 4-7/8" x 4-7/8" x 6-1/8"



50Z459

**Latex-Free Insti-Dam®
20pk**
\$69

Ideal for anterior isolation. Made with polyisoprene.

Dimensions: 4-1/4" diameter



Lingua-Fix®

Large suction surface avoids blockage and reduces sedimentation of drilling byproducts in pipes and tubing.

- Single-Use Only
- Comes pre-assembled

50Z953S

**Lingua-Fix®
10pk**
\$16

50Z953

**Lingua-Fix®
50pk**
\$57



comfortable,
hands-free suction

Pink Petal™

- Holds low-volume suction tip in place
- Protects and retracts the cheek
- Patient can close around suction with ease
- Easily slides onto saliva ejector
- Comfortable material
- Single-Use Only
- Latex-free

50Z976

Pink Petal™ 100pk
\$54



CRYSTALHD®

What a view!

The brighter the workspace, the better the illumination and visibility. With Zirc's Crystal HD® lens, you are visually working in a brighter environment.


If you haven't made the change to HD – what are you waiting for?

As a microscope user, I can understand the importance of clear vision. Zirc mirrors significantly improved the brightness of my image, giving my eyes highest quality vision. It's true – brightness makes all the difference.

Dr. Dimitris Tsanaktsidis
Microscopic & Aesthetic Dentistry


Why Use Zirc's Resin Mirrors?

Crystal HD® reflects color 99% accurately while Rhodium only reflects color 75% accurately.





Brighter than rhodium mirrors!

Image provided by Dr. Lori Trost



	Resin	Stainless
Autoclavable	X	X
Lens flush with mirror frame	X	
No galvanic shock on patient	X	
Rust-free	X	

CUTTING EDGE IMAGERY



What's Your Mirror Size?

diagram is to scale

#0 Mirror
14mm#3 Mirror
20mm#4 Mirror
22mm#5 Mirror
24mm

ZIRC TIP

Crystal HD® Mirror Heads**Crystal HD® Resin Mouth Mirror
Screw-on Heads - Cone Socket 12pk
\$85****Single Color**Available in all
colors and K (black)#4 | 50Z371
#5 | 50Z372**Assorted Vibrant**2 of each color:
[N, O, P, Q, R, S]#4 | 50Z371-NEON
#5 | 50Z372-NEON**Assorted Jewel**2 of each color:
[E, J, T, I, M, U]#4 | 50Z371-JWL
#5 | 50Z372-JWL**Assorted Classic**2 of each color:
[A, B, D, G, I, K]#4 | 50Z371-PSTL
#5 | 50Z372-PSTL**Crystal HD® Mouth Mirrors****Crystal HD® Mouth Mirror
Soft Grip 12pk
\$85****Assorted Vibrant**
2 of each color
[N, O, P, Q, R, S]#4 | 50Z363-NEON
#5 | 50Z364-NEON**Assorted Jewel**
2 of each color
[E, J, T, I, M, U]#4 | 50Z363-JWL
#5 | 50Z364-JWL**Crystal HD® Mouth Mirror
Soft Grip Black/Gray**

\$32 4pk

#4 | 50Z361
#5 | 50Z362

\$85 12pk

#4 | 50Z363
#5 | 50Z364**Crystal HD® Mouth Mirror
Thin Grip 12pk
\$85****Assorted Vibrant**
2 of each color
[N, O, P, Q, R, S]#4 | 50Z353-NEON
#5 | 50Z354-NEON**Assorted Jewel**
2 of each color
[E, J, T, I, M, U]#4 | 50Z353-JWL
#5 | 50Z354-JWL**Crystal HD® Resin Mouth Mirror
Double-Sided Heads - Cone Socket 4pk
\$46****Assorted Vibrant**1 of each color
[N, O, R, S]

#4 | 50Z374-NEON

Gray Mirror Head

#4 | 50Z374

**Crystal HD® Stainless Mouth Mirror
Screw-On Head - Cone Socket
\$85****Stainless Steel 12pk**#3 | 50R379
#4 | 50R380
#5 | 50R381**Stainless Steel Head**

\$14 Each

#0 | 50R388



Autoclavable



Washer Disinfect



Surface Disinfect

Mirror Handles



50Z342

Ergo Grip Handle - Cone socket

\$18.50 Each Available in K (black) & all colors except A



50Z340

Soft Grip Handle - Cone Socket

\$18.50 Each Available in all colors



MEGA Mirrors

* NOT a medical device

50Z309

MEGA Mirror

\$22 Color code: N, P, Q, R, S

A fun toy/prop for your patients.
2-1/2" lens and 15-1/2" in length.

S

Mirror Magic® Anti-Fog System



50R201

Mirror Magic® Organizer

\$16.50
Holds: 48 Swipe Pads and 1-6ml bottle of Anti-Fog Solution. Store on operator counter or in drawer. Non-slip bottom for secure placement.

Dimensions: 3-3/4" x 3-5/8" x 1-7/8"

Note: box only, does not include Mirror Magic®

S



50R205

Mirror Magic® Anti-Fog Solution 10pk

\$76
Each bottle contains 6ml each and lasts approx. 50 patients.



50R210

Mirror Magic® Swipe pads 96pk


\$27 Single-Use Only



50R200

Mirror Magic® Trial Kit

\$26
Comes with 2-6ml bottles of Anti-Fog Solution and 24 Swipe Pads.



Denture & Retainer Boxes



ZIRC TIP



- 1 Effortlessly upload your logo or choose from pre-made icons to imprint with our high quality, 4-color printing process.
- 2 Personalize with up to 4 lines of text
- 3 Manufactured and imprinted in the USA!

Retainer & Denture Boxes S

One piece construction. Solid locking mechanism.

Retainer Box 12pk

1" Deep | 25R500
\$14

Dimensions: 3-1/8" x 2-3/4" x 1"

1 1/2" Deep | 25R550
\$15

Dimensions: 3-1/8" x 2-3/4" x 1-1/2"

Single Colors



Assorted Pack



30R800
Denture Box 12pk
\$16

Dimensions: 4-3/8" x 3-3/8" x 2"

Single -or- Assorted Pack



Personalized Imprinting

Upload your own logo or select a full color icon from our pre-made fun or professional sets. Add up to four lines of text for your business name, phone number, website, or your own unique message. It's as easy as that!

Imprinted Retainer Box

24pk

1" Deep | 25P500
\$40

Dimensions: 3-1/8" x 2-3/4" x 1"

1 1/2" Deep | 25P550
\$42

Dimensions: 3-1/8" x 2-3/4" x 1-1/2"

1000pk

1" Deep | 25PB500
\$1,351

Dimensions: 3-1/8" x 2-3/4" x 1"

1 1/2" Deep | 25PB550
\$1,422

Dimensions: 3-1/8" x 2-3/4" x 1-1/2"

Single Colors



Assorted Pack



Imprinted Denture Box

24PK

30P800
\$44

Dimensions: 4-3/8" x 3-3/8" x 2"

1000PK

30PB800
\$1,460

Dimensions: 4-3/8" x 3-3/8" x 2"

Single Colors



Assorted Pack





The Steri-Soaker S

It is designed to hold instruments for easy cleaning and is large enough to accommodate Zirc Steri-Containers, Cassettes, Bur Blocks, and Bur Guards.

50Z805
STERI-SOAKER

\$107 Color code: B, I

Calendar on lid conveniently reminds staff when to change the solution; use color code rings to mark dates. Transfer basket's ribbed and perforated design allows quick drain-off and easy pick-up of instruments. Holds 3.78L of liquid; fill lines shown in quarts and liters. Hinged fitted cover.

Dimensions: 13-7/8" x 9-5/8" x 5-1/8"

Safe to use chemicals: cleaners (enzymes, detergents), disinfectants (phenols)/cold sterilants for glutaraldehyde and chlorines.



50Z417
CAD/CAM Block Organizer

\$55 Color Code: E, J, N, Q, S, T, U

Provides a professional platform to present and discuss same-day restoration. Fits all CAD/CAM milling blocks including PrograMill®, CEREC®, PlanMill™ Blocks. Store in drawer or cabinet.

Dimensions: 5-9/16" W x 1-7/8" x 7-3/4" D

S

Instrument Mats



20Z480
Mini Matt
\$15.50 Available in all colors


Holds 8 instruments. Organizes instruments at chairside and for sterilization. Non-slip feet for stability.

Dimensions: 3-7/8" x 3-1/2" x 1/2"

50Z917
Instrument Strap
\$6.50/Pair Gray only

Locking strap secures instruments.


Dimensions: 3-5/8" x 1/2" x 5/16"



20Z464
Instrument Matt
\$17 Available in all colors

Holds 10 hand instruments. Reverse side holds 4 surgical instruments. Non-slip, soft material. Organizes instruments in order of procedure use. Can be cut to fit Mini Tray.

Dimensions: 6-3/4" x 2" x 3/8"



E-Z Storage S

E-Z Storage Tub System

Hinges allow cover to stay open on its own! Each tub comes with one vertical divider and one half-depth adjustment divider. Hinged cover comes in white or clear.

29R605 E-Z Storage Organizer (White Cover)

\$38 Rack sold separately

Dimensions: 5-7/8" x 4-7/8" x 5-3/4"

29R610 E-Z Storage Organizer (Clear Cover)

\$38 Rack sold separately

Dimensions: 5-7/8" x 4-7/8" x 5-3/4"

29R627 Replacement Hinges (Pair)

\$5.50

Dimensions: 7/8" x 1/2" x 1/8"

Hinges allow cover
to stay open



Rack holding angled Tub

29R650 E-Z Storage Single Rack **\$26**

Holds storage organizer Tub. Tub can be placed flat or at an angle for easy access of supplies. Wall or cabinet mounted with optional mounting bracket. Sturdy metal construction; plastic coated. White, non-slip feet.

Dimensions: 7" x 5" x 4-3/4"

29R651 E-Z Storage Mounting Bracket **\$5.50**

White metal bracket for the E-Z Storage Tub single rack. Pre-drilled holes for cabinet or wall mount.

Dimensions: 5-1/8" x 1" x 1/2"



Handi-Hopper

Handi-Hopper (With Mounting Bracket) S

Small plastic waste receptacle for disposal of contaminated waste at chairside, keeping tray clean. Slotted mounting bracket (included) for easy removal and replacement. To use for infectious waste, order red (Biohazard) liners. To use for other waste, order clear liners.

20Z404 Handi-Hopper

\$22 Color code: G, I

Dimensions: 7" x 4" x 6-3/8"

20Z408 Handi-Hopper Bracket

\$4.50 Color code: G, I

Dimensions: 2-1/4" x 1/4" x 1"

Mount this mini waste receptacle within arm's-reach of chairside for convenient disposal of contaminated waste.

Handi-Hopper Liners **\$25**

20Z405 | 100/BOX
Red Biohazard

Single-Use Only

20Z410 | 100/BOX
Clear

Single-Use Only



22Z100 Waste-Drop Ring **\$16** Color code: A, G, I

Waste-Drop Ring for cabinet top. Fits 5" hole. Made of self-lubricating, easy-to-clean resin.

Dimensions: 5-5/8" Dia. x 3"

S

E-Z Access Shelf & Accessories



Provides convenient work space for auxiliary items

20Z420

E-Z Access Shelf

\$130 Color code: A, G, I

E-Z mounting to standard equipment post (1-3/4" x 2-1/4"); No tools needed. Will not damage post's finish. E-Z adjustment to multiple positions. Contoured edge to contain spills and materials. Stain resistant finish; Surface disinfect only. Weight capacity of 25lbs. Rotates 360°.

Standard E-Z Access Shelf as is accommodates 1-3/4" x 2-1/4" post.

Workable tray area: 13-3/4" L x 10-1/4" W



E-Z Access Shelf Adapters & Mounts S



20Z422

E-Z Access Shelf Adapter for Smaller Posts

\$7 Color code: A, G, I

Adapts to posts 1-3/8" to 2-3/4"

Dimensions: 2-1/8" x 1-7/8" x 3/4"

20Z437

E-Z Access Shelf 3 1/2" Bolt Adapter Kit

\$12 Color code: A, G, I

Accommodates post sizes from 2-1/4" to 2-3/8"

Dimensions: 4" x 2" x 2"

20Z438

E-Z Access Shelf 4 1/2" Bold Adapter Kit

\$12 Color code: A, G, I

Accommodates post sizes from 2-3/8" to 2-3/4"

Dimensions: 4-1/2" x 2" x 2"



20Z436I

E-Z Access Shelf Cup Holder Snap-On

\$12 Gray only

Fits all standard cups. No tools needed.

Dimensions: 3-1/8" x 2-1/2" x 1"



20Z428

E-Z Access Shelf 8" Extension Arm

\$35 Color code: A, G, I

Enables shelf to extend over patient. Shelf has strength to hold 25lbs. when connected to arm. May be used with Wall Mount, Surface Mount, or directly to post.

Dimensions: 10" x 4" x 4-1/2"



20Z426

E-Z Access Shelf Wall or Cabinet Mount

\$32 Color code: A, G, I

Mount on vertical flat surface only.

Dimensions: 6-1/2" x 5-3/4" x 4-7/8"




20Z424

E-Z Access Shelf Surface Mount

\$17 Color code: A, G, I

Mount on horizontal flat surface only.

Dimensions: 3-3/4" x 3-3/4" x 2-7/8"



Water Wise S

Water Wise vs. Traditional Water Bottle

Features	Water Wise	Water Bottle
Easy to clean	X	
Antimicrobial	X	
Easy-to-fill large mouth	X	
Quick release	X	
Level indicator coded for: Sterile and distilled water, Medicament	X	
Strong	X	X
Easy to grip	X	
Sealable lid	X	
No threads to strip	X	
Pressure relief valve	X	



23Z100
Water Wise Starter Kit
\$36

Holds 1 liter. Includes bottle, adapter, lid, neon yellow and neon blue bubbles (water level indicator), and tube.

23Z105
Water Wise Replacement Bottle
\$24

Holds 1 liter. Includes bottle, blue lid, and neon yellow and neon blue bubbles (water level indicator), and tube.

- 1 Quick release adapter Pressure relief valve
- 2 Water level indicator (bubbles)
- 3 Flat stable bottom with large opening for easy cleaning and filling

Dispensers S



28R800
Cotton Roll Dispenser
\$34 Color code: A, B

Heavy-duty, no-skid base. Dispenses up to three rolls at a time. Extra sturdy construction.

Dimensions: 4-3/16" x 2-1/16" x 4"



55Z200
Matrix Dispenser
\$28.50 Color code: A, B

One roll of 3/8" x 12 yards Matrix clear tape 0.002 gauge. Allows matrix to be contoured (u-shaped) as it dispenses. Stainless steel cutter allows for desired length. Disposable.

Dimensions: 2-1/8" x 2-1/4" x 5/8"



Pellet Dispenser
\$34 White only

Spring loaded and pushes pellets to the surface.

28R825
Pellet Size 3, 4, and 5
Dimensions: 1-7/8" x 1-7/8" x 2-1/16"

28R850
Pellet Size 0, 1, and 2
Dimensions: 1-7/8" x 1-7/8" x 2-1/16"



Mighty Mixer Bowls, Spatulas, and Organizers

Mighty Mixer Spatulas

Self-lubricating for easy cleaning. Broad, medium flex blade. Ideal for alginates and plastic.



50Z501
Single Spatula
\$6.50 Each White only

Dimensions: 7-5/16" x 3/16" x 1-1/4"



50Z501ASTD
Assorted 4 Pack
\$22

One (1) of each color: J, N, P, and R

Dimensions: 7-5/16" x 3/16" x 1-1/4"



50Z523
Delrin Spatula 10pk
\$11

Flexible flat blade. Non-stick. Resistant to chemicals. Great for mixing cements, acrylics, etc.

Dimensions: 4-5/8" x 3/16" Dia.



50Z502
Cement Spatula
\$7 Each

White, plastic, double-end spatula. For fillings, cement, and pastes. Self-lubricating plastic for easy cleaning.

Dimensions: 6-3/8" x 1/4" x 3/8"



50Z521
Mixing Sticks 100pk
\$24 White only

Single-Use Only

Dimensions: 4-1/2" long



20Z496V
Mixing Well 144pk
\$26 White only

Single-Use Only

Dimensions: 4-1/2" long



Mighty Mixer Bowls 36pk
\$26

Disposable, eliminating cleanup time. Fits securely into mixing bowl without slippage. Saves wear on mixing bowl. Free-standing for use with or without mixing bowl. Capable for use with electric Alginate mixers.

Single-Use Only

50Z505 | White
50Z505ASTD | Assorted

Dimensions: 4-7/8" Dia. x 3-1/8"



20Z490B
Implant Organizer
\$29 Blue only

Quick and easy recognition. Individual compartments for implants. Shape of the organizer matches the arch of the mouth for implant placement. Separate compartments prevent implant appliances from mixing together. Lubricated plastic enables tooth number to be written in the designated area. Cover helps keep materials secure.

Dimensions: 6-7/8" x 4-5/8" x 5/8"



20Z495R
Restorative Organizer
\$34 Neon purple only

Compartments shaped like the arch of the mouth for veneers, inlays, composite, and onlays. Disposable Mixing Wells included for bonding agent and cement. Accessory area for brushes and materials. Lubricated plastic enables tooth number to be written in the designated area. Light protective cover can be used as a light shield and to prevent materials from setting early. Includes 25pk of Mixing Wells.

Dimensions: 6-7/8" x 4-5/8" x 5/8"





Every office can improve on efficiency.





We can help.

Our clinically-trained consultants are ready to help you lead the way in transforming your practice into the most efficient version of itself.





ZIRC.COM | (763) 682-6636

*All U.S. pricing subject to change